



SAN LUIS OBISPO COUNTY Public Health Laboratory

Phone: (805) 781-5507 Fax: (805) 781-1023 www.sloPublicHealth.org/lab

LITERATURE SUMMARY

TB NAAT performance

The table summarizes the sensitivity and specificity of the Amplified Mycobacterium tuberculosis Direct Test (MTD, Hologic-Gen-Probe) and the Gene Xpert MTB-RIF assay (Cepheid, Inc) for detection of *M. tuberculosis* nucleic acid in primarily respiratory specimens.

NAAT	Specimens n=	Sensitivity	Specificity	Ref
MTD	1,151 resp.	91.7% overall 98.7 % smear pos 62.2 % smear neg	98.7 % overall 97.8% smear pos 98.9% smear neg	1
MTD	162 resp	97.3% overall	87.1.% overall	2
Gene Xpert	+non-resp	90.9% overall	89.0%	2
Gene Xpert	107 resp	100% smear pos 71.7% smear neg	100%	3
Gene Xpert	85 resp	75.3% smear neg	100%	4
MTD	31 resp 4 resp	93.5% smear pos 75% smear neg		5
Gene Xpert	217 resp	98% smear pos 72% smear neg		6
Gene Xpert	27 35	100% smear-pos 68.6% smear neg		7

References:

1. R L. Guerra, MD; N M. Hooper, JF. Baker, et al Use of the Amplified *Mycobacterium tuberculosis* Direct Test in a PublicHealth Laboratory-Test Performance and Impact on Clinical Care
Chest 132: 946-951. 2007.
- 2 J Teo, R Jureen, D Chiang, D Chan and R Lin. Comparison of Two Nucleic Acid Amplification Assays, the Xpert MTB/RIF Assay and the Amplified Mycobacterium tuberculosis Direct Assay in for Detection of Mycobacterium tuberculosis in Respiratory and Non-respiratory Specimens.
J. Clin. Microbiol., 49:3659. 2011.
3. D Helb, M Jones, E Story, et al. Rapid Detection of *Mycobacterium tuberculosis* and Rifampin Resistance by Use of On-Demand, Near-Patient Technology. J. Clin. Microbiol., 48:229. 2010.
4. R Moure, L Munoz, M Torres et al. Rapid Detection of *Mycobacterium tuberculosis* Complex and Rifampin Resistance in Smear-Negative Clinical Samples by Use of an Integrated Real-Time PCR Method. J. Clin. Microbiol., 49:1137. 2011.

LITERATURE SUMMARY

TB NAAT performance

5. L Han, PEIvin, and JBernardo. Nonclinical Selection Criteria for Maximizing Yield of Nucleic Acid Amplification Tests in Tuberculosis Diagnosis. *J. Clin. Microbiol.*, 50:2592. 2012.
6. E Marloew, S Novak-Weekley, J Cumpio, et al. Evaluation of the Cepheid Xpert MTB/RIF Assay for Direct Detection of *Mycobacterium tuberculosis* Complex in Respiratory Specimens. *J. Clin. Microbiol.*, 49:1611. 2011.
7. A Zeka, S Tasbakab, and C Cavusoglu. Evaluation of the GeneXpert MTB/RIF Assay for Rapid Diagnosis of Tuberculosis and Detection of Rifampin Resistance in Pulmonary and Extrapulmonary Specimens *J. Clin. Microbiol.*, 49:4138. 2011.